



صنایع پزشکی

سرما بهبود

سپنتا

SBS Medical Ind.



دانش بنیان

Knowledge Enterprise



MANA مانا
Neonate Cooling
سردکننده نوزادان

Model : **SBNeo**

Total body cooling device for treatment of neonates suffering hypoxic ischemic encephalopathy (HIE)

The first and sole
licensed manufacturer in Iran

Sarma Behboud Sepanta Medical Industries
Mashhad Health Science and Technology Park,
Bu-Ali Square, Mashhad, Iran

+98 51-35015023

0915 098 9047

www.sbmed.ir

info@sbmed.ir

sbmed.ir

sbmed.ir



IMED
National Medical
Device Directorate IR, IRAN



Treatment Mechanism




By placing the neonate on a cold mattress, the core temperature decreases to 33.5 degrees Celsius

The core temperature is maintained within this range for 72 hours

Decrease in temperature slows down the metabolic rate

Brain cells are able to recover, which ultimately reduces or stops further brain damage

Total Body Cooling Timeline



Functional Features

■ Two main Modes of Operation include:

a) Single-point temperature control mode

In this mode, the clinician can set the rectal target temperature and the device will automatically perform cooling or heating operations at the tuned rate or at maximum power, depending on user's choice. Then the device will maintain the target temperature without interruption and with minimum fluctuations. This operation mode effectively transforms the device into a temperature management system that can be used for hypothermia, hyperthermia, or normothermia treatments.

b) Multi-point temperature control mode

This operation mode is based on the standard protocol of asphyxia treatment. The user determines the target rectal temperature (by default 33.5 degrees Celsius), and the device cools the infant in the shortest possible time. The core temperature is then maintained during the requested period, which is 72 hours in the standard protocol. Finally, the device gradually rewarms the neonate at a set rate (maximally, 0.5°C/hour) until reaching the normal temperature of 37 degrees or any other desired temperature.

- Simultaneous reading and display of surface and core body temperature using skin and rectal sensors
- The temperature display accuracy is 0.1°C, the maximum temperature stabilization error is 0.5°C, and temperature sampling is done every 1 second
- The possibility of viewing changes in rectal temperature during the entire treatment period as a graph with various graphic capabilities
- Easy and safe use of machine due to design of the device's accessories storage space inside the
- Advanced and adaptive temperature control algorithm with self-tuning capability
- The possibility of immediate fast rewarming during treatment at the discretion of the treatment staff
- Multilayer hardware and software temperature protection systems

● Technical Specifications

Dimensions

50*40* 98 CM

Weight **52 Kg**

Coolant Temp

16-42 °C

Noise intensity

53db at 1 meter

Filling Jar volume

0.5 Lit

Power supply

220 Volt 50 Hz

Power consume

520 Watt

Coolant tank volume

3.5 Lit



Sensors

- Reusable paediatric rectal probe
- Disposable paediatric rectal probe
- Reusable paediatric skin probe
- Disposable paediatric skin probe

Connecting pipes

- Silicone hoses with flex insulation
- Protected with fabric cover

Cooling System

- Thermoelectric modules

Structure

- Metal structure with electrostatic heat treated paint

User interface

- 7-inch touch LCD

Heating system

- Electric element heaters

Connectors

- Non-spill quick couplings

Controller

- PID Dual mode

Coolant

- Pure water

Mattress

- Biocompatible
- Washable
- Uniform heat treatment

Approvals and Certificates

